Universe of **Biocoral**[®]


Calcium
Carbonate
naturel
wholly
mineral

*“Full fadom five thy father lies
of his bones are coral made
those are pearls that were his eyes
nothing of him that doth fade
but doth su er a sea-change
into something rich and strange.”*

William Shakespeare
THE TEMPEST (Act 1, SC II).

Biocoral® is the only natural wholly mineral bone graft substitute having a highly porous architecture, which is formed from naturally occurring aragonite crystal with three-dimensional interconnectivity allowing for optimum formation of new bone.

Biocoral®:

- Natural Calcium Carbonate (CaCO_3),
- 3 dimensional randomly interconnected pore,
- Pore size: $150\mu\text{m}$,
- Porosity: 50% allowing an ideal newly bone formation.

THE PHYSIOLOGY OF BONE GRAFTING

The biology of bone grafts and their substitutes is appreciated from an understanding of the bone formation processes of Osteogenesis, Osteoinduction and Osteoconduction.

Osteogenesis: The cellular elements within a donor graft, which survive transplantation and synthesize new bone at the recipient site.

Osteoinduction: New bone is formed through the active recruitment of host mesenchymal stem cells from the surrounding tissue, which differentiate into bone-forming osteoblasts. This process is facilitated by the presence of growth factors within the graft, principally Bone Morphogenetic Proteins (BMP).

Osteoconduction: The facilitation of blood-vessel incursion and new-bone formation into a defined passive trellis structure.

All bone graft and bone-graft-substitute materials can be described through these processes.

CHARACTERISTICS OF BIOCORAL[®]

Biocoral[®] is a natural bone substitute used since more than 30 years in all surgical, repair and bone regeneration procedures. Biocoral[®] is available in the form of granules, beads, blocks, and the shaped prostheses. The Biocoral[®] is the only natural wholly mineral resorbable bone substitute composed of Calcium Carbonate (>98%).

Biocoral[®] biocompatibility, together with its osteoconductive and osteophilic properties, induces specific biological activity in the recipient bone, similar to the physiological natural bone metabolism. This activity leads to graduate resorption of Biocoral[®] by osteoclasts and its replacement by osteoblasts in newly formed bone.

The porosity architecture of Biocoral[®] is in the form of aragonite crystalline. The regular structures of pores, their volume, their size and the thickness of their walls are the special characteristics of Biocoral[®]. The open porosity of Biocoral[®], allows the blood cells and bone marrow cells (blood, anions, cations, etc...) to spread and infiltrate in its core which speed up the bone ingrowths. Biocoral[®] present the remarkable mechanical resistance qualities associated to its porosity (50% and 20%) similar to those of the cancellous and cortical bone respectively.

Section of a human femoral cortex

Section of Biocoral[®]

GOOD REASONS TO USE BIOCORAL[®]

Biocompatible

Biocoral[®] is perfectly tolerated by the human body with no risk of contamination and is compatible with the structural requirements of bone growth.

Bioresorbable

Biocoral[®] has excellent bone integration with total resorption between 3 to 9 months.

Biocoral[®] because of its mineral and architectural characteristics (Aragonite Crystal and Porosity), once placed in bony sites is quickly impregnated with autogenous blood or bone marrow, with a proved calcification from the day 9th.

Osteoconductor

Biocoral[®] has an ideal porosity which allows a quick invasion of bone marrow and integration of newly formed bone.

Re-Initiates Bone Mineralization Process

Biocoral[®] is used as an active ingredient for reinitiating the process of bone remineralization.

Replaced by Newly Formed Bone

Biocoral[®] is quickly vascularized and progressively resorbed by osteoclast cells, which is then replaced by osteoblast cells in order to conduct newly formed bone identical to the recipient bone.

No Risk of Viral Transfer & Contamination

Biocoral[®] follows strict quality control procedures which are performed at each stage of manufacturing process and guarantees its compliance with high quality standards and offers surgeons a truly trustworthy biomaterial.

Easy to Use

Biocoral[®] is available in a variety of shapes and sizes. It can easily be used and applied to the surgical site once infiltrated by autogenous blood or bone marrow.

Biocoral[®] is available in different forms: granules, beads, blocks and shaped prostheses.


Avoids taking autologous bone graft & Reduce Cost

Biocoral[®]'s uses avoids taking autologous bone graft (which can lead to risks such as postoperative infection, bleeding, pain) and prevent unwanted aesthetic disorders at donor site. The take of autologous bone graft increases the time of the surgical process and the costs related into. (the surgical act and miscellaneous costs).

Biocoral[®] is the best alternative to the autologous grafts and also contributes to decrease and minimize patient's hospitalization and its health care costs.

CHEMICAL COMPOSITION OF BIOCORAL®

Following the studies of Pr Le Petitcorps in 2006 at Hospital University Center of Bordeaux, L'ICMCB-ENSCP (one of the CNRS laboratories (UPR 9048)), Biocoral®'s chemical composition is confirmed being wholly mineral as described below:

The image shows a table titled "Chemical composition of Biocoral®" with a background of a microscopic view of bone tissue. The table lists various components and their percentages, separated by horizontal lines.

<i>Chemical composition of</i> Biocoral®	
Calcium Carbonate (CaCo3)	> 98 %
(including calcium)	> 40 %
Trace elements (including Fluorine and Strontium)	0.7 à 1 %
Magnesium	0.05 à 0.2 %
Sodium	< 1 %
Potassium	< 0.03 %
Phosphore (phosphate form)	< 0.05 %
Water	< 0.5 %

Several components are present at levels equivalent to those found in mammalian bone, notably trace elements, which play a vital role in the process of mineralization and in the activation of enzymatic reactions in bone cells.

Two trace elements have specific effects:

- Strontium is involved in the formation and growth of the crystalline component of bone, protects calcification mechanisms and increases mineralization. Furthermore, Strontium levels are higher in the most active bone structure the metaphysis and bone callus.
- Fluorine in proper quantity increases bone formation by direct effects on proliferation of cellular precursors of osteoblasts.


TECHNICAL PROOF ON BIOCORAL[®]

A. Determination of Pore Size


Biocoral[®] pore size is important as this offers the ideal environment for vascularization and migration of osteoclasts and osteoblasts. It ranges between 20% and 50% depending on the species selected. Acropora (20% of porosity) is very close to the cortical bone. The Porites (50% of porosity) is very close to the cancellous bone. That porosity allows the surgeon to choose Biocoral[®], according to the clinical indication. However, those specifications do not constitute a rigid standard. Natural Calcium Carbonates with different porosities can be used according to the operating procedure.

Biocoral[®] pore size ranges between 150 to 500 microns, depending on the species, selected according to clinical indications. It has previously been shown that these sizes are optimal for occupation by fluids and bone marrow cells in order to complete mineralized newly formed bone. For indications where dense material is necessary, for example for strong mechanical compression strains, the natural microporous Calcium Carbonate (Acropora) can be used.

Bone cells (bone marrow and blood of the recipient bone) can freely invade the open porous structure of Biocoral[®] deep in its core. This cellular invasion determines the first phase of the bone restoration process characterized by the development of a neo-vascularization. Some corals (Porites in particular) have the architecture similar with the cancellous bone.

Section of Biocoral[®] 50% porosity

Section of Biocoral[®] 20% porosity

B. Overall porosity



Biocoral[®] has a maximum porosity of 50% for allowing the ideal bone ingrowths. To induce the bone regeneration process overall porosity is one of the important key factors and its high porosity enhances the osteoconductivity.

MECHANICAL CHARACTERISTICS	BIOCORAL [®]		
	DENSE	20% POROSITY	50% POROSITY
Breaking deformation	0.47 ± 0.03 %	0.53 ± 0.12 %	0.26 ± 0.04 %
Breaking stress in compression in MPa	395 ± 29 330 % that of fresh cortical bone 50 % that of titanium alloy	110 ± 32 75 % that of fresh cortical bone	25.2 ± 5.2 25 % that of fresh cortical bone 500 % that of cancellous bone
Young's modulus (modulus of elasticity) in GPa	101.6 ± 3.7	24.6 ± 3.30	7.99 ± 0.37

(Orthopedics Researches Laboratory, U.A. C.N.R.S., 1161,
 Expert evaluation report, INOTEB, 1987)

C. Three-Dimensional interconnected pore for ideal Bony Growth

Biocoral[®]'s architecture is entirely porous and is defined by the total volume three-dimensional randomly interconnected to allow bone formation throughout the entire implant. (Photos below)

BIOCORAL[®]'S INDICATIONS

Biocoral[®] is a natural bone substitute used since more than 30 years in all surgical, repair and bone regeneration procedures.

Biocoral[®] is the only biocompatible and bioresorbable calcium salt used as an active ingredient for local treatment of diseases associated with demineralization or mineralization defects of bone, with the aim of reinitiating the bone remineralization process (Patented Application).

Orthopedic surgery:

- Substitution of cancellous bone graft,
- Insufficiency of autogenous bone graft,
- Filler of metaphyseal bone graft (osteoporotic disease),
- Femoral neck fracture associated with a cancellous bone demineralization,
- Demineralized bone cavity filling (bone-cyst, chondroma, acetabular cavity),
- Upper and lower limb fractures associated with bone demineralization,
- Non-union atrophic bone fractures,

Cranio-Maxillo-Facial surgery:

- Maxillo reconstructive and plastic surgery,
- Plastic and facial reconstructive surgery of Malar,
- Nasal sulcus filling,
- Plastic and reconstructive surgery of Orbit (post-traumatic, post-cancerous),

Oral surgery:

- Filling, reconstruction of bone defects and regeneration of bone for: Periodontal, implant and endodontic surgery,
- Filling, reconstruction of bone defects and regeneration of bone for: Pre prosthetic surgery and implant procedure: Filling after tooth extraction,
- Oral surgery: Filling of cystics cavities.


Particular Applications

- Sinus lift augmentation,
- Lateral sinus lift augmentation,
- Filling of voluminous cystics cavities,
- Filling and augmentation of lateral sinus.

Biocoral[®] is a bone graft substitute which provides an osteoconductive matrix and it is a very useful for adding volume to autogenous graft to obtain osteoinduction. It is recommended to impregnate Biocoral[®] with bone marrow or blood, particularly when it is necessary to obtain immediate cohesion between Biocoral[®] granules or beads. Once blood is mixed with Biocoral[®], the existence of fibrin in blood makes a composite easier to fill in the recipient bone.

HISTOLOGIC ANALYSIS OF FUSION MASS

A. From Biocoral[®] to Newly Formed Bone

Phase 1: Invasion by blood cells and extravasated bone marrow.


Phase 2: Vascularization (v: vessels).

Phase 3: Resorption of Biocoral[®] by osteoclasts (oc: osteoclast).


Phase 4: Bone neoformation with osteoblastic (ob: osteoblast) apposition and concomitant resorption.

Phase 5: Remodeling of neoformed tissue to produce the architecture of the recipient bone.

B. New bone formation associated with the bone's physiological modification

BIO : Biocoral[®] NB : Newly Formed Bone
 → : Osteoblast cells


CT : Connective Tissue
 → : Osteoclast cells


We can notice the physiological bony modification with an osteoblastic phase and an osteoclastic phase.

RADIOLOGICAL RESULTS OF FEMORAL LENGTHENING


An eight year old child had a right leg that was more than five centimeters shorter than the left leg. Such a discrepancy normally calls for an elongation of the malformed leg. A Biocoral[®] graft was implanted. During the following years, Biocoral[®] was perfectly integrated into the site and there was a successful reconstruction of the femoral shaft.

Post Operative

One Year Later

Three Years Later

REFERENCES BIBLIOGRAPHICS

- ✓ **FUNDAMENTAL RESEARCH**
- ✓ **CLINICAL RESEARCH - ORTHOPAEDICS AND SPINE SURGERY**
- ✓ **CLINICAL RESEARCH - MAXILLOFACIAL SURGERY**
- ✓ **CLINICAL RESEARCH - ENT SURGERY**
- ✓ **CLINICAL RESEARCH - ORAL SURGERY**

REFERENCES BIBLIOGRAPHICS

FUNDAMENTAL RESEARCH

PETITE H.P., VIATEAU V., POLLACK C, BENSARD W, MEUNIER A, BOURGUIGNON M, OUDINA K, SEDEL L, GUILLEMIN G. (LRO-CNRS, Université Diderot, Paris, FRANCE)

Tissue-engineered bone regeneration
Nature Biotechnology, 2000; 18: 959-963

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Evaluation of combinations of Titanium, Zirconia, and Alumina Implants with 2 Bone Fillers in the Dog.
Int. J. Oral&Maxillofacial Implants, 1999;14:271-277

GAO T.J, TUOMINEN T.K., LINDHOLM T.S., KOMMONEN B., LINDHOLM T.C.

(University of Tampere and Helsinki, FINLAND)

Morphological and biomechanical difference in healing in segmental tibial defects implanted with Biocoral[®] or tricalcium phosphate cylinders.
Biomaterials, Elsevier Science Limited., 1997; 18: 219-223,

FUJIMORI Y., SUGAYA K., IZUMI H., KOZAWA Y. (Nihon University, Matsudo, Chiba, JAPAN)

Scanning Electron Microscopic observation of Bone Marrow-derived Osteoclast-like cells resorbing Coral.
J. Oral Biol., 1997; 39: 241-246

REIS S.A., VOIGT C., MULLER-MAI C HERBST H, BISSON S, GROSS U. (Free University of Berlin, GERMANY)

Procollagen $\alpha 1$ (I) transcripts in cells near the interface of coralline implants in rats, detected by in situ hybridization
Clin Oral Impl Res, 1996; 7: 253-260

MOON I.S., CHAI J.K., CHO K.S., WIKESJÖ U., KIM C.K.

(Yonsei University, Seoul, KOREA)

Effects of polyglactin mesh combined with resorbable calcium carbonate or replamineform hydroxyapatite on periodontal repair in dogs.
J. Clin. Periodontal, 1996; 23: 945-951

BRAYE F., IRIGARAY J.L., JALLOT E., OUDADESSE H., WEBER G., DESCHAMPS N., DESCHAMPS C., FRAYSSINET P., TOURENNE P.,

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Diffusion of mineral elements evaluated by PIXE at the bone-coral interface.
J. Biomater. Sci Polymer Edn, 1996; 7 (n°8): 741 - 749.

MULLER-MAI C., VOIGT C., DE ALMEIDA REIS SR, HERBST H, GROSS U.M. (Free University of Berlin, GERMANY)

Substitution of natural coral by cortical bone and bone marrow in the rat femur. Part II Sem, Tem and in situ hybridisation.
Journal of Materials Science : Materials in Medicine, 1996; 7: 479 - 488.

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Resorption of natural calcium carbonate by avian osteoclasts in vitro.
Cells and Materials, Vol 5, N° 2, Pages 157 - 165, 1995.

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Comparison of the ossification kinetics after implantation of a radioactivated coral and a natural coral.
Journal of Materials Science : Materials in Medicine 6 (1995) 230-234.

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Engineering bone: challenges and obstacles
J. Cell Mol. Med. Vol. 9, No 1, 2005 pp. 72-84

KANIA RE, MEUNIER A, HAMADOUCHE M, SEDEL L, PETITE H.
(Laboratoire de Recherches Orthopédiques, Université D. Diderot, Paris VII, URA CNRS 1432, Paris, France)
Addition of Fibrin Sealant to Ceramic Promotes Bone Repair : Long-Term Study in Rabbit Femoral Defect Model
Journal of Biomedical Materials Research. 1998 Spring;43(1):38-45.

PETITE H.P., BLANCHAT C., TRIFFITT J.T, OUDINA K, SEDEL L, BENSALD W. (UMR. -CNRS, Faculté de médecine Lariboisière saint – Louis, université D. Diderot). Biomaterials 11 December 2002
A biodegradable fibrin scaffold for mesenchymal stem cell transplantation

PETITE H.P., BLANCHAT C., OUDINA K, SEDEL L, BENSALD W, D.V.S.,M.S.,VIATEAU V., Ph.D. , POTIER E., BOUSSON V., M.D., B.S., GUILLEMIN G.
De Novo Reconstruction of functional bone by tissue engineering in the metatarsal sheep model
The journal of Tissue Engineering .vol. 11, Number 5/6, 2005

VIATEAU V, GUILLEMIN G, BOUSSON V, OUDINA K, HANNOUCHE D, SEDEL L, LOGEART-AVRAMAGLOU D, PETITE H.(Ecole nationale vétérinaire d'Alfort, Laboratoire de recherches orthopédiques, Université D. Diderot, laboratoire de radiologie expérimentale)
Long- bone critical- size defects treated with Tissue- engineered Grafts: A study on sheep
Journal of Orthopedic Research society 2007, published by Wiley InterScience,

CLINICAL RESEARCH - ORTHOPAEDICS AND SPINE SURGERY

E. SOFFER¹, J.P. OUHAYOUN², A. MEUNIER¹, F. ANAGNOSTOU²
Effects of Autologous Platelet Lysates on Ceramic Particle Resorption and New Bone Formation in Critical Size Defects:
The Role of Anatomical Sites

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Received 25 March 2005; revised 8 September 2005; accepted 26 October 2005.

Published online 16 March 2006 in Wiley InterScience (www.interscience.wiley.com). DOI: 10.1002/jbm.30516

W. BENSALD^{a,b}, J.T. TRIFFITT^b, C. BLANCHAT^a, K. OUDINA^a, L. SEDEL^a, H. PETITE^{a,*}

A biodegradable fibrin scaffold for mesenchymal stem cell transplantation

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^b Botnar Research Laboratory, Nuffield Department of Orthopaedic Surgery, University of Oxford, UK

Received 5 September 2002; accepted 11 December 2002

W. BENSALD, K. OUDINA, V. VIATEAU, E. POTIER, V. BOUSSON, C. BLANCHAT, L. SEDEL, . GUILLEMIN, and H. PETITE, De Novo Reconstruction of Functional Bone by Tissue Engineering in the Metatarsal Sheep Model

TISSUE ENGINEERING, Volume 11, Number 5/6, 2005

© Mary Ann Liebert, Inc.

D. LOGEART-AVRAMAGLOU, F. ANAGNOSTOU, R. BIZIOQ, H. PETITE*

Engineering bone: challenges and obstacles

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Received: December 22, 2004; Accepted: January 25, 2005

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Tissue-engineered bone regeneration

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Osteogenesis with coral is increased by BMP and BMC in a rat cranioplasty

^a Unite& de chirurgie craniofaciale, Département de neurochirurgie, Hopital Necker, 75743 Paris cedex 15, France

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^c Sulzer Orthopedics, 4056 Youngxeld Street, Wheat Ridge, CO 80033, USA

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A physiological approach in stabilization and consolidation of unstable femoral neck fracture in osteoporotic elderly patients: a retrospective review
Eur J Orthop Surg Traumatol; 2003; 13:145-155

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Behavior of natural coral in a human osteoporotic bone

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ZAIOUR W., DEHOUX E., DEPREY F., SEGAL Ph.
(*Centre Hospitalier Universitaire de Reims, FRANCE*)
Use of coral as a bone graft substitute for anterior fusion of the lower cervical spine. A review of twenty cases.
Orthopaedic Product News, May/June 92

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KOROVESSIS.P, KOUREAS.G, ZACHARATOS.S, PAPAZISIS.Z, LAMBIRIS.E
Correlative radiological, self-assessment and clinical analysis of evolution in instrumented dorsal and lateral fusion for degenerative lumbar spine disease. Autograft versus coralline hydroxyapatite] *Eur Spine J* (2005) 14 : 630-638

RAMZI.N, RIBEIRO-VAZ.G, FOMEKONG.E, LECOUVET.F, RAFTOPOULOS.C.
Long term outcome of anterior cervical discectomy and fusion using coral graft] *Acta neurochir* (2008) 150 : 1249 – 1256

NICOLAIDES AP, PAPANIKOLAOU A, POLYZOIDES AJ.

Successful treatment of valgus deformity of the knee with an open supracondylar osteotomy using a coral wedge: a brief report of two cases. *Knee*.2000 Apr 1;7(2):105-107

ROUX.FX, BRASNU.D, MENARD.M, DEVAUX.B, NOHRA.G, LOTY.B.
Madreporic coral for cranial base reconstruction]. *Acta neurochir (Wien)* (1995) 133 : 201-205

CIROTTEAU Y.
The behavior of natural coral in a context of diaphyseal atrophic pseudoarthrosis. *Eur J Orthop Surg Traumatol*, 2004, 14 : 89-98. France

CLINICAL RESEARCH - MAXILLOFACIAL SURGERY

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SOUYRIS F., PELLEQUER C., PAYROT C., SERVERA C.
(*University Hospital Centre Montpellier, FRANCE*)
Coral, a New Biomedical Material, Experimental and First Clinical Investigations on Madreporaria.
max.fac. Surg. 1985; 13: 64-69

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PETITE H.P., SEDEL L, MEUNIER A., DAMIEN C., DE POLLAK, ARNAUD E.
(*LRO, Université D. Diderot*), *Biomaterials* 16 May 1999
Osteogenesis with coral is increased by BMP and BMC in a rat cranioplasty

BAE Y.C., CHOI S.J., MOON J.S., NAM S.B.
Comparison of the postoperative outcome in pure medial orbital fracture among three groups: using porous polyethylene or hydroxyapatite through subciliary approach and transnasal endoscopic correction.
Ann Plast Surg. 2007 Sep;59(3):287-90.

NAM S.B., BAE Y.C., MOON J.S., KANG Y.S.
Analysis of the postoperative outcome in 405 cases of orbital fracture using 2 synthetic orbital implants. *Ann Plast Surg.* 2006 Mar; 56(3):263-7

DAGLI A.S., AKALIN Y., BILGILI H., SECKIN S.
Correction of saddle nose deformities by coral implantation. *Eur Arch Otorhinolaryngol.* 1997;254(6):274-6

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CARLO MANGANO. (Department of Biomaterials Science, Università dell'Insubria, Varese, Italy)
Combining Scaffolds and Osteogenic Cells in Regenerative Bone Surgery : A Preliminary Histological Report in Human Maxillary Sinus Augmentation. Carlo Mangano, DDS, MD; Alessandro Mangano, MD; Giovanna Iezzi, DDS, Fabio L. Borges, DDS ; Susana d'Avila, DDS, MS, PhD; Jamil Awad Shibli, DDS, MS, PhD.
Clinical Implant Dentistry and Related Research, Volume, Number, 2009

YUKNA R. A. (Louisiana State University, USA)
A 5 year follow-up of 16 patients treated with coralline calcium carbonate (Biocoral™) bone replacement grafts in infrabony defects.
J. Clin Periodontol 1998; 25: 1036-1040

CORRENTE G. (University of Turin, ITALIA)
Supracrestal Bone Regeneration around dental implants using a Calcium Carbonate and a Fibrin-Fibronectin sealing system: Clinical and Histologic evidence. *Int. J. Periodontics & Restorative Dentistry*, 1997; 17, 171-181.

OUHAYOUN J.P, (Université Paris VII-Garancière, FRANCE)
Bone grafts and Biomaterials used as bone graft substitutes
II European Workshop in Periodontics, Switzerland, Quintessenz Verlag, 1996; 313-348

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MORA F, OUHAYOUN J.P (Université Paris VII-Garancière, FRANCE)
Clinical evaluation of natural coral and porous hydroxyapatite implants in periodontal bone lesions:
Results of 1 year follow-up.
J. Clinical Periodontol 1995; 22: 877-884

YUKNA R. A. (Louisiana State University, USA)
Clinical evaluation of coralline calcium carbonate (BIOCORAL) as a bone replacement graft material in human periodontal osseous defects.
Journal Periodontol, February 1994, Vol 65 - Nb. 2, pp 177-185.

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
Biocoral®

Indications

Biocoral® is used as bone graft substitute in all surgical, repair and bone regeneration procedures.

Biocoral® is the only biocompatible and bioresorbable calcium salt used as an active ingredient for local treatment of diseases associated with demineralization or mineralization defects of bone, with the aim of reinitiating the process of bone remineralization. (Patented Application)

It is appropriate to read and follow the recommendations described in the leaflet accompanying the medical device Biocoral®.

Advantages

- Natural wholly mineral biomaterial,
- Easy to use and easy handling,
- No risk of viral transfer & contamination,
- Quick infiltration by autogenous blood or bone marrow once placed in bony sites,
- Variety of shapes and sizes for many fields.

Biocoral® has remarkable physical, chemical and architectural properties similar to those of the human bone.

General Data

Biocoral® is a medical device classified in the class III according to the European regulation. Biocoral® is in complete conformities with the European and International norms. Biocoral® although marketed since the end of the years 80, was the first bone substitute registered in France with TIPS (Tarifs Interministériels des Prestation Sanitaires). Nevertheless within the framework of its registration, Biocoral® has been received approving opinion of Microbiological Committee of Health on July 05, 1995 referenced under the number 9600201B01.

Biocoral® obtained the authorization of marketing in the European countries market "EC lable" on December 30, 1996.

The manufacture of Biocoral® comply with the regulation inforc.

Biocoral®'s radiosterilization for pharmaceutical use is obtained by β . Rays with the delivery dose between 25 and 50 KGrays. The shelf life of Biocoral® is five years from the date of sterilization.